


Mobilett Plus/Plus E/Plus M

SP



Installation

Installation Instructions

Mounting of giraffe stickers

© Siemens AG 1998

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Chapter	Page	Revision
All	All	01

	Page
1 Prerequisites	1 - 1
General	1 - 1
Tools required	1 - 1
Time required	1 - 1
Parts required	1 - 1
2 Location of stickers	2 - 1
3 Installation procedure	3 - 1
Installation of sticker 1-7 (See fig 2.1)	3 - 1
Installation of sticker 8-11 (See fig 2.1)	3 - 2
Final procedure	3 - 2

This page intentionally left blank.

General

The application surfaces should be clean and dry before the assembly.

The stickers should not be exposed to abnormal stretching or pulling. This could result in bubbles or a deterioration of the shape stability.

Avoid touching the adhesive side of the sticker. This could lead to a decreased adhesive capacity of the sticker.

The best application temperature is 15-20°C, but you can still apply the stickers at temperatures as low as 10°C.

If air-bubbles will appear you can puncture these with a needle and release the air with help of the plastic applicator or with the finger.

To keep the stickers as fresh as possible it is important to keep in mind that the stickers have their optimal adhesive period at approximately 48 hours after the assembly. Because of this you should wait at least a couple of days before you wash the stickers the first time.

Tools required

- Masking tape

Time required


Approximately one hour for 1 person.

Parts required

Installation kit, part no. 64 58 942, containing the following parts:

Qty	Part
1	Installation Instructions
1	Giraffe toy
1	Plastic applicator
11	Stickers

This page intentionally left blank.



MOB00033

Fig. 2.1 The numbers of the stickers and the placing of the tape

The installation kit of giraffe stickers contains 11 different stickers. See fig. 2-1. The figure also shows where to place the masking tape.




Figure 2.2 shows where the stickers numbered 1-6 are to be pasted on the MOBILETT Plus / Plus E / Plus M.





Figure 2.3 shows where the stickers numbered 7-11 are to be pasted.


Installation of sticker 1-7 (See fig 2.1)

NOTICE


While attaching sticker 8 make sure that the red focus indicator is visible.



First you attach the stickers with masking tape to the covers as shown in fig 2.2 and 2.3. The tape will be used as a hinge. In figure 2.1 you can see where to place the tape.



Lift the sticker all the way to the tape-hinge at an angle of 180°. Remove the release paper and carefully lower the sticker, while you paste the sticker to the cover with the other hand or with the plastic applicator.




Use the plastic applicator from the hinge to the side. Work with overlapping moves. Remove the tape and follow the same procedure on the other side of the sticker. If air-bubbles will appear you can puncture these with a needle and release the air with help of the plastic applicator or with the finger.

Installation of sticker 8-11 (See fig 2.1)

Remove the release paper and carefully paste the sticker to the covers according to fig 2.3. If air-bubbles will occur follow the procedure as described above.

Final procedure



Finally remove the protective paper at an angle of 180°.